

7"

3"

**Use:** SoSweet is a convenient, flavored syrup based vehicle formulated to assist the compounding pharmacist in the preparation of oral extemporaneous delivery systems. SoSweet may be used alone, to flavor aqueous solutions, or in conjunction with suspending agents to impart flavor and sweetness while increasing the physical stability of the suspension. SoSweet contains suitable preservatives and is buffered to a slightly acidic pH to provide a vehicle which is dependable, consistent pharmaceutically elegant and easily utilized by the compounding pharmacist.

**Directions:** SoSweet may be used alone or as a flavored syrup vehicle. SoSweet may be diluted to desired taste or thickness with water or other aqueous solutions. SoSweet is also designed to complement the use of any suitable Suspending Vehicle. As such, it is flavored and buffered to be added to a pre-formed suspension and diluted in approximately 50-50 concentrations.

**FOR BEST RESULTS:** Aqueous Solutions: Wet the material with water or other water-miscible substances such as glycerin or sorbitol. Alternately, dissolve the material in minimal amounts of water, add SoSweet to dilute to volume. If a suitable suspending agent is used an approximate equal portion of each is optimal.

NDC #31722-959-01

**SoSweet** 

**Flavored  
Syrup Vehicle**

16 fl. oz. (1PT.) 473 ml

**Contains:** Purified water, Sucrose, Glycerin, Sorbitol and flavoring. Buffered with Citric acid and Sodium phosphate. Preserved with methylparaben and potassium sorbate.


**PRECAUTION:** SoSweet is formulated to a slightly acidic pH to enhance the stability of compounds subject to degradation at higher pH's. Basic, ionized drugs with low solubility in the unionized state may precipitate. If such a compound is encountered or anticipated, a suspension might be considered for these compounds. As with all extemporaneous formulations, no guarantees regarding the suitability or stability of your unique formulation can be made. Without specific information regarding your compound, "Keep Refrigerated" and "Shake Well" labeling and conservative expiration dating should be employed.

Distributed by:


**CAMBER**  
PHARMACEUTICALS, INC.

Piscataway, NJ 08854



-  Black Print
-  Reflex Blue
-  PMS# 200 Red
-  Die  
(Does Not Print)

**#4 out/left**



### IMPORTANT NOTICE

**The approval of all label proofs is requested within TWO business days of receipt to ensure a quick turn time. After this time delays may occur with production.**

I have carefully read the proof and authorize you to proceed as indicated.

Signed: \_\_\_\_\_ Date: \_\_\_\_\_

☐ 1ST REVISION

☐ 2ND REVISION

☐ 3RD REVISION

- ☐ OK'D As Is
- ☐ OK'D With Corrections
- ☐ Request Another Proof